

Алматы (7273)495-231	Иваново (4932)77-34-06	Магнитогорск (3519)55-03-13	Ростов-на-Дону (863)308-18-15	Тольятти (8482)63-91-07
Ангарск (3955)60-70-56	Ижевск (3412)26-03-58	Москва (495)268-04-70	Рязань (4912)46-61-64	Томск (3822)98-41-53
Архангельск (8182)63-90-72	Иркутск (395)279-98-46	Мурманск (8152)59-64-93	Самара (846)206-03-16	Тула (4872)33-79-87
Астрахань (8512)99-46-04	Казань (843)206-01-48	Набережные Челны (8552)20-53-41	Санкт-Петербург (812)309-46-40	Тюмень (3452)66-21-18
Барнаул (3852)73-04-60	Калининград (4012)72-03-81	Нижний Новгород (831)429-08-12	Саратов (845)249-38-78	Ульяновск (8422)24-23-59
Белгород (4722)40-23-64	Калуга (4842)92-23-67	Новокузнецк (3843)20-46-81	Севастополь (8692)22-31-93	Улан-Удэ (3012)59-97-51
Благовещенск (4162)22-76-07	Кемерово (3842)65-04-62	Новобрянск (3496)41-32-12	Саранск (8342)22-96-24	Уфа (347)229-48-12
Брянск (4832)59-03-52	Киров (8332)68-02-04	Новосибирск (383)227-86-73	Симферополь (3652)67-13-56	Хабаровск (4212)92-98-04
Владивосток (423)249-28-31	Коломна (4966)23-41-49	Омск (3812)21-46-40	Смоленск (4812)29-41-54	Чебоксары (8352)28-53-07
Владикавказ (8672)28-90-48	Кострома (4942)77-07-48	Орел (4862)44-53-42	Сочи (862)225-72-31	Челябинск (351)202-03-61
Владимир (4922)49-43-18	Краснодар (861)203-40-90	Оренбург (3532)37-68-04	Ставрополь (8652)20-65-13	Череповец (8202)49-02-64
Волгоград (844)278-03-48	Красноярск (391)204-63-61	Пенза (8412)22-31-16	Сургут (3462)77-98-35	Чита (3022)38-34-83
Вологда (8172)26-41-59	Курск (4712)77-13-04	Петрозаводск (8142)55-98-37	Сыктывкар (8212)25-95-17	Якутск (4112)23-90-97
Воронеж (473)204-51-73	Курган (3522)50-90-47	Псков (8112)59-10-37	Тамбов (4752)50-40-97	Ярославль (4852)69-52-93
Екатеринбург (343)384-55-89	Липецк (4742)52-20-81	Пермь (342)205-81-47	Тверь (4822)63-31-35	

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

<https://bandelin.nt-rt.ru/> || bni@nt-rt.ru

SONOREX TECHNOLOGY MB 75.2 LIFTING DEVICE



The lifting device makes work easier for the operator. It enables automatic oscillation into the cleaning or rinsing bath and facilitates the transport of the filled product baskets from tub to tub. The slow oscillating up and down movement increases the cleaning effect of the ultrasound and rinses off the dissolved dirt more effectively.

-The filled cleaning basket is hung manually in the lifting device. By switching the timer, the lifting device lowers the basket into the tub and automatically starts oscillation. After the timer has expired, the oscillation stops and the lifting device automatically lifts the basket out of the tub back into the starting position.

ORDER NUMBER

8309

OPERATION

Display

TIME SETTING [MIN]

1-99, infinite

OPERATING VOLTAGE [V]

400 V

MAINS CABLE

5

OUTER DIMENSION LENGTH [MM]

depending on the number of tubs

OUTER DIMENSION WIDTH [MM]

740

OUTER DIMENSION HEIGHT [MM]

1700

DEGREE OF PROTECTION HOUSING

IP 20

PAYLOAD [KG]

40

PACKING LENGTH [MM]

0,8

PACKING WIDTH [MM]

0,6

PACKING DEPTH [MM]

0